





harex

Medical Device for management of male urinary incontinence

www.harex.net







https://www.facebook.com/harexdevice https://www.twitter.com_harexdevice https://www.instagram.com/harexdevice

Distributor




Manufacturer



harex

Medical Device for management of male urinary incontinence



harex

is the solution to improve your quality of life.





Male urinary incontinence is recognized by the World Health Organization as a chronic disease, which decisively affects people's quality of life in social situations, as well as physically, sexually, psychologically and professionally.

HAREX is a medical device that is placed just below the head of the penis and by gentle pressure on the urethra, stops the flow of urine, preventing and controlling involuntary leakage of urine.

It's small, comfortable and fully adjustable!

Comfort Harex is designed to be worn for many hours without discomfort. You won't notice it, and they won't notice it.

> **Durable and reliable** Manufactured with high quality materials that guarantee perfect functioning over a long period of time.

> > Easy to use You can place and adjust the device without anyone helping you, so you can be completely autonomous.

% plastic materials Made entirely of sanitary latex-free plastic material.

Easy cleaning Just apply warm water and soap. It does not contain metallic pieces which can oxidize or generate dermatological problems.

Allows blood circulation The points needed to apply pressure to the urethra have been calculated, so that only the urine outlet duct is blocked, thus allowing free blood circulation.



Solution

More than 400 million men worldwide suffer from urinary incontinence. **Harex** is the solution that can help them return to a completely normal life.



Quality of life

As it is very easy to use and intuitive, this means you can regain your day-to-day freedom, and improve your quality of life.



Trust

Using **Harex**, you will be able to carry out any sport, work or leisure activity with physical effort, without the inconvenience. You will never go through embarrassing situations again.




Recovery of normality

Enjoy the daily activities you gave up again; walk, run, travel, work, ... enjoy total normality with your partner, family and friends.



Easy to put and use.

- down. Loosen the belt.
- **2**-Introduce the penis into the Device and adjust the belt to the
- **3** Move the pressure activator up, to fix the belt.



urethra in such a way as

Belt Manufactured

discomfort.

LATEX

FREE

with lightweight and high quality sanitary materials,

which exert the correct

pressure without causing

to prevent involuntary urine leakage.

Belt adjustment lever

Precise adjustment

Ridged belt for convenient

positioning and adjustment.

on the penis.

A safety system that allows the belt to be fixed in the optimum position and ensures control by means of a fitted closure.

