

CODE	DESCRIPTION	SPECIMEN/SAMPLE SITE
M451101	TEMPLATE - Sampling Template - 10x10 cm / 100 sqcm - Individually Wrapped	Environmental control sampling template

INTENDED USE

The template is designed to be used as a sampling guide for all surface samples that require a 10x10cm sample area.

INSTRUCTIONS FOR USE

- Open the rip seal on the plastic bag.
- Holding the template handle with the plastic bag. (You can also use a sterile glove if preferred)
- Fold the bag back on itself over your hand.
- Place on the surface to be sampled.
- Discard template after use

PRECAUTIONS

Single use, disposable swabbing template.

STORAGE

Keep dry and store at ambient temperature


EXPIRY DATE

5 Years shelf life

SYMBOLS APPLICABLE TO FAMILY GROUP



SYMBOLS & DEFINITIONS



Temperature Limitation



Use By (YYYY/MM)



Method of Sterilisation Using Irradiation



Manufacturer



Do Not Reuse



Consult Instructions for Use



LRQA Approval Mark



Date of Manufacture



Do Not Use If Package Damaged



Peel Here



Contains Sufficient for <n> Tests



Reference /Catalogue
Number



LOT or batch identification