

Pathfinder

Used for continuous flow and pulsatile action with minimal tubing, ensuring efficient and controlled irrigation.



Key Features

Improved Visualization

- ▷ Ensures a clear view by quickly clearing debris from the surgical field
- ▷ Laser beam scatter is reduced, providing accurate delivery of energy

Selective Irrigation

- ▷ Slow continuous flush controllable by saline bag height and/or pressure cuff
- ▷ Gentle squeeze clears large amounts of debris from the field

Control

- ▷ Pathfinder plus bulb fits in scope hand, preserving dominant hand for instrumentation

Universally Compatible

- ▷ Standard Luer-lock fittings attach directly to any Cystoscope or Endoscope

Manufactured by:

Marflow AG, Alte Landstrasse 54, CH-8546. Islikon, Switzerland
Tel. +41 44 709 01 01 | Fax. +41 44 710 21 36
www.marflow.ch
marflow@marflow.ch